LC121/c - 06/2015

The Bluetooth® word mark and logos are registered trademarks owned by Bluetooth SIG, Inc. and any use of such marks by Laica S.p.A. is under licence.
Other trademarks and trade name are those of their respective owners.
iOS, iPhone, iPad mini, iPod touch, App Store are trademarks of Apple Inc., registered in the U.S. and other countries - Android™, Play Store™ and Google Play™ are trademarks of Google Inc.

LAICA S.p.A.
Viale del Lavoro, 10 - Fraz. Ponte
36021 Barbarano Vicentino (VI) - Italy
Tel. +39 0444.795314 - 795321
Fax. +39 0444.795324
Made in China

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GB Instructions and guarantee

SMART ELECTRONIC SCALE WITH BODY FAT %, WATER %,
MUSCLE %, BONE MASS AND BODY MASS (BMI) INDICATOR

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This new electronic scale lets you calculate your body composition (percentage of body fat, water and muscle and BMI calculation) based on your height, age and gender. Exact knowledge of your body composition (calculation of lean mass, fat mass, total body water and muscle mass) enables a correct assessment of excess fat and at the same time helps identify the desirable weight based on nutritional condition. In this way, a proper diet can be adopted along with physical activity as fundamental conditions for improving quality of life, promoting health and preventing disease. Capacity 180 kg / 397 lb / 28 st 5 lb - division 100 g / 0.2 lb.

- Attention!**
- **The body fat, water and muscle percentages along with bone mass and BMI should only be considered to be indicative. For further information, consult a physician or dietician.**
 - **Self-measurements are for monitoring purposes, not diagnosis or treatment. Unusual values must always be discussed with your physician. Under no circumstances should you change the dosage of any medications prescribed by your physician.**

BEFORE USING, READ ALL OF THE INSTRUCTIONS AND WARNINGS CONTAINED IN THIS MANUAL AND KEEP IT IN A SAFE PLACE.

BODY COMPOSITION
It is important to know individual body composition in order to be able to evaluate health. In fact, weight is not enough data. The quality of the weight must be defined, which means dividing the human body into two large masses as follows:

- lean mass: bone, muscle, viscera, water, organs, blood;
- fat mass: fat cells.

At equal weights, these two masses take up different volumes. Through the electrodes located on the platform, the scale sends a weak, low intensity electrical current that passes through the body. The electric current reacts differently depending on the mass it encounters. Lean mass provides little resistance while fat mass provides greater resistance. This body resistance, called bioelectrical impedance, also varies depending upon gender, age and height. For proper monitoring of your ideal weight, it is important to check not only weight but also to compare the body fat percentage calculate by the scale with the parameters given in the table and asses the proper percentage of water and bone mass with the information given below.

If the percentages are outside of normal values, we recommend consulting a

physician. We recommend keeping track of these values for a certain period of time in order to be able to understand them better.

For a valid comparison, you should always weigh yourself at the same time of day and under the same conditions. The best time to weigh yourself is in the morning before eating, when undressed. Wait about 15 minutes after rising so water can distribute itself through your body.

body fat %
The body fat percentage indicated in the table below represents the total fat that is given by adding essential (or primary) fat and storage fat. Essential fat is indispensable as it is necessary for the metabolism. The percentage of this fat differs for the two genders: About 4% for men and about 12% for women (% calculated on total body weight). Storage fat, which is in theory superfluous, is necessary in minimum quantities for good health, especially for those who perform physical activities. In this case, the percentage of this fat also differs for the two genders: 12% men, 15% women (% calculated on total body weight). So, the normal body fat percentage of body weight is 16% for men and 27% for women. These percentages vary with the age of the individual as shown in the table below.

BODYFAT TABLE					
	Age	% Underweight Low	% Normal Normal	% Slightly overweight High	% Overweight V.High
WOMEN	<39	<21	21.2-33	33.1-39	>39
	40-59	<23	23.1-34	34.1-40	>40
	>59	<24	24.1-36	36.1-42	>42
MEN	<39	<8	8.1-20	20.1-25	>25
	40-59	<11	11.1-22	22.1-28	>28
	>59	<13	13.1-25	25.1-30	>30

muscle %
Muscles are the skeleton's motors and by means of the energy released during contraction, they allow us to move, manage and control the movements of the head, trunk and upper and lower limbs. There are three different types of muscle: smooth, skeletal and cardiac. Cardiac muscle is found in the heart. This is the only organ made of striated fibres and that contracts involuntarily.

Instructions and guarantee

Smooth muscle is made of smooth fibres. These are found exclusively in the walls of the digestive tract, respiratory tract and genital tract and in the wall of the venous and arterial blood vessels. Smooth muscle is under involuntary control. Skeletal muscle, which is called voluntary because it contracts according to our will, includes the muscles attached to the skeleton and the facial muscles. These voluntary striated muscles move the skeleton, thus allowing us to move and communicate with the environment outside our bodies. This scale indicates the % of skeletal muscle and smooth muscle tissue in the body. On average, 40% of body mass is made up of skeletal muscle. In general, the percentage can increase in athletes during prolonged physical activity. In moderate and daily physical activity, the % tends to increase, but slowly (the difference is seen after months of activity). Physical activity also improves the water balance in the body. The percentage of muscle starts to decrease with age (after 40). Men usually have more accentuated muscles than women.

Water %
Water is the main component of lean mass. The percentage of water in lean mass is 73%. Lean mass is obtained by subtracting an individuals body fat mass from the total weight (example: 25 year old man with fat mass 17% = lean mass 83%). About 60% of body weight is water, which is the basic constituent of all living things and is essential to life. Women (who normally have more body fat then men) have a water percentage of around 55%-58% while men have around 60%-62%. The percentage of water is greater in childhood (at birth it varies from 70%-77%), decreases with age (in the elderly the percentage varies from 45%-55%) and as fat deposits increase. Obese individuals have less water mass than lean individuals. Refer to the data below to compare the percentage detected. Water performs many biological functions:

- it promotes digestion;
- it is a solvent for substances that are important for the body;
- it removes metabolic waste;
- it transports oxygen, nutrients, enzymes, hormones and glucose to the cells;
- it is a source of mineral salts;
- it regulates body temperature.

On a daily basis, the human body expels a certain amount of water which must be regularly replaced. A loss of water of about 10% of body weight can compromise health.

Body weight variations can occur in a short period of time (hours), a medium period of time (days) or over a long period (weeks). Small variations in weight over a few hours is normal and due to water retention since extracellular water (which includes the interstitial fluid, plasma, lymph and transcellular fluid) is the only compartment capable of varying so rapidly. When weight decreases in a few days, it is probably a decrease in total water (about 60% inside the cells - intracellular water - and about 40% outside the cells - extracellular water) and thus even of lean mass. Body fat changes are ultimately subject to variations in weeks. **Therefore, when assessing your body composition trend, it is necessary to keep the following in mind:**

- **A decrease in weight accompanied by no change in body fat or even an increase in body fat means that your body only lost water (for example, after exercising, time in the sauna or a limited diet for rapid weight loss).**
- **An increase in weight accompanied by a decrease or no change in body fat indicates the development of precious muscle mass.**
- **When the weight and body fat decrease at the same time it means that your diet is working and volume is being lost.**
- **A high percentage of body fat or reduced muscle mass always indicates a body water percentage that is too low.**

Warning!
When following a weight loss treatment, body fat should decrease and not muscle mass.

The skeletal system
The skeleton is the sturdy and flexible structure that supports our body, giving it shape, protecting the internal organs, anchoring the muscles and allowing movement. It is primarily made up of bones (206) and cartilage. Most of the bones in the female skeleton are the same as those in the male skeleton, except smaller. The exception is the pelvis. In fact, the female pelvis is usually larger than the male pelvis. This scale gives the weight of the skeleton. The skeletal weight reading must be related to the age, height and lifestyle of each individual.

Body mass index
The body mass index (BMI) is the relation of body weight to height. The classification given be low is correlated to the risk that body weight exerts on health. When BMI is between 18.5 and 24.9 it means that the weight does not statistically impact

health and is to be considered normal. Any deviation, above or below this BMI interval should be considered a health risk. In this case, we recommend consulting a physician. Getting in the habit of following a balanced diet definitely helps improve quality of life, promote health and prevent illness, increase vitality and maintain ideal weight!

BMI TABLE	BMI VALUE
UNDERWEIGHT	<18.5
NORMAL	18.5-24.9
OVERWEIGHT	25.0-29.9
OBESE	>30

Examples of body composition

Person 1
Data entered: female, height 163 cm, age 37.
Data calculated by the scale: Weight 52.9 kg, Body fat 21.2%, Body water 57%, Muscle 35.2%, skeletal weight 4.7 kg, BMI 19.9.
Body composition analysis:
The table shows that this is a normal body composition.
Let's take an in-depth look.
The lean mass percentage is 78.8% (100%-21.2%).
The percentage of water is 57% (78.8% x 73%).
The body fat is as follows:
• essential fat => 52.9 kg x 12%= 6.3 kg
• storage fat => 52.9 kg x 15%= 7.9 kg
The total body fat (essential fat + storage fat) is given from 6.3+7.9 kg = 14.2 kg which is 27% of the total body weight (52.9 kg) (ideal body fat for person 1).
The difference between 21.2% (% real body fat %) and 27% (ideal body fat %) is -5.8%, which when applied to the total weight of 52.9 kg is equal to 3.06 kg.
3.06 kg is the number of kg that person 1 could increase to obtain ideal body fat.
The water percentage at 57% is within the standard that women should have.
The percentage of the skeletal and smooth muscles (35.2%) is within the standard (women usually have less accentuated muscles than men).

- TIPS FOR USE
- Place the scale on a flat, sturdy surface.

An uneven or broken floor (e.g., rug, carpet, linoleum) could affect the accuracy of the weight.

- Step on the scale with your feet parallel and your weight equally distributed (see drawing). Wait without moving until the scale computes your weight.
- It is important to weigh yourself with dry, bare feet positioned correctly on the electrodes (see figure). Also avoid contact between the legs so the scale operates correctly.
- The scale is not dangers since it operates with disposable, 1.5 V batteries. Nevertheless, it is not recommended for use by the following people:
 - individuals wearing a pacemaker or other electronic medical device
 - people with fever
 - children under 10 years of age
 - pregnant women
 - individuals undergoing dialysis, with oedema in the legs, dimorphism or osteoporosis
 - people who take cardiovascular medications
- Always weigh yourself on the same scale, every day at the same time, preferably undressed and before eating breakfast. To optimize the scale, weight yourself twice and if the two weights are different, your weight is between the two. Wait about 15 minutes after rising so water can distribute itself through your body.
- The surface of the scale can be slippery when wet.
- If there is an error while analysing your body fat, water and muscle percentage, the scale will automatically shut off.
- Scale reading precision could vary in the following cases:
 - high level of alcohol
 - high level of caffeine or drugs
 - after rigorous sports activities
 - during illness
 - during pregnancy
 - after heavy meals
 - if you have dehydration problems
 - if your data (height, age, gender) were not entered correctly
 - if your feet are wet or not clean
 - people with swollen legs

Attention!!
The body fat, water and muscle percentages along with bone mass and BMI should only be considered to be indicative.
For further information, consult a physician or dietician.

Instructions and guarantee

INSTRUCTIONS FOR USE

a. Inserting batteries
This scale operates with 4 disposable 1.5V AAA alkaline batteries. When the batteries are dead, the battery symbol will appear on the display. Remove the batteries if not using the scale for long periods of time. Remove dead batteries by removing the lid from the battery compartment on the back of the device. Remove the dead batteries and dispose of them as special waste at those locations indicated for recycling. For more information regarding the disposal of used batteries, contact the store where you purchased the appliance that contained the batteries or the City or local waste disposal service.
Attention: For greater safety, removal of the batteries by children under 12 years of age is NOT recommended.

b. "UNIT" key
This key is located on the bottom of the scale and lets you choose the unit of measure to use when weighing: kg / lb / stone. When the scale is off or on, select the unit of measure to use (kg, lb or st) with the "UNIT" button on the bottom of the scale.

c. Operation as a scale only
1) With the scale off, step on with your feet parallel and your weight equally distributed (see drawing).
Attention: only turn on after stepping on the scale. Wait without moving until the scale computes your weight. Do not lean on anything during weighing to avoid altering the reading.
2) Step off the scale. The display will continue to show the weight for a few seconds and then will shut off automatically. Already starting from the second use of each memory area, the scale can recognize the memory area of each of the 8 users by their weight (for weight changes within 2 kg).
In this case, when the scale shows your weight, it automatically displays the calculation of the % of body fat, % of water, % of skeletal musculature and weight of the skeletal system (see paragraphs "Enter your personal data" and "Calculation of weight, % of body fat, % of water, % of skeletal musculature and weight of the skeletal system").
When the scale does not recognize your memory area, it suggests in any case the memory areas where similar data already exist.
At this point, press the keys "◀" or "▶" to select your memory area (fig. 1).

d. Entering personal data
This scale can memorize data for 8 different people.
1) With the scale off, press the "SET" button (on the side above the display). The scale will emit a short beep.
2) The display will show "P1" which indicates person 1. If you wish to change the memory area (from 1 to 8), press the "◀" or "▶" buttons. Once you have chosen the memory area, confirm by pressing "SET". The scale will beep each time the button is pressed.
3) The display will now show the shape of a man, woman, a running man (🏃) and a running woman (🏃). Use the "◀" and "▶" keys to choose and confirm by pressing SET.
4) Now enter your height (from 100 to 220 cm) using the "◀" or "▶" symbols to increase or decrease the value. The display will give the value in cm or ft. Press "SET" to confirm.
5) Use the "◀" and "▶" symbols again to enter your age (from 10 to 100) and confirm by pressing "SET".
6) The scale has now saved all of your data. At this point the scale shows "0.0" and will turn off after a few seconds.
To change the data you entered, recall the memory area and repeat the entire procedure. To insert data for another person, choose a free memory area and repeat the entire procedure given above.
Attention: when you are in the "personal data insertion" function the display will show "0.0" if no button is pressed within 10 seconds. The scale shuts off automatically after a few seconds.

e. Calculate weight, body fat %, water % and muscle %
1) With the scale off, press the "SET" button.
2) The display shows "P1".
At this point press the "◀" or "▶" button until the desired memory area is shown. After a few seconds the display will show "0.0".
3) Step on the scale, making sure your bare feet are correctly positioned above the two electrodes. The display will show the following data in sequence 3 times:

- weight
- the percentage of fat in your body (👤) and the comparison with the parameters in the "body fat table" (Low, Normal, High, V.High)"
- the percentage of water (💧)
- the percentage of muscle (💪)
- the bone mass (🦴)
- the "📱" symbol: this symbol will disappear as soon as the scale sends your data to the smartphone. If the symbol remains on the display, it means that the data were not sent to the smartphone. They remain saved in the scale until sent again to the smartphone

(see paragraph "Operation with app").
When calculating body composition, the scale shows "□□□□".
Compare the body fat percentage with the parameters given in the "body fat table". If it is outside normal values, consult your physician. Also assess the water percentage as shown in the paragraph "water %" and the muscle percentage as in the paragraph "muscle %".
4) The scale will shut off automatically after the third time the data above has been displayed.
Attention! If the display does not show "0.0" it will not be possible to read your weight and calculate the fat and water percentage. If over the allowed weight, the display will show "----". In this case, remove the weight right away to avoid damaging the scale mechanism.

- f. Operation with app**
List of compatible devices:
- **iOS:**
 - iPhone 6 Plus
 - iPhone 6
 - iPhone 5s
 - iPhone 5c
 - iPhone 5
 - iPhone 4s
 - iPad (4th generation)
 - iPad (3rd generation)
 - iPad mini
 - iPod touch (5th generation)
 - For Android™ devices, check the bluetooth version, downloading the free app BLE CHECKER from Play Store™.

Support bluetooth 4.0 technology, you can send your body composition data (fat, water and muscle percentage and BMI) to your smartphone to better analyse your weight trend.
During the following operations, we recommend keeping the smartphone near the scale to keep the connection active.
1) Activate bluetooth in the smartphone settings.
2) Install the free app "LAICA" in your smartphone.
3) Start the app and follow the instructions to create a new user, entering your personal data.
4) To connect the scale to your smartphone, click on the icon in the upper left on your phone (👤) and select "Device". Now click on "Add device" and hold down the "UNIT" button on

- the bottom of the scale until the scale display shows "0 0" and select the scale when it is recognized by your phone.
- 5) The app will now ask you to select one of 8 available memory areas. This way, the scale will only send the data related to the selected memory area to your smart phone ("P1", "P2", "P3", etc.).
- 6) Now you can view and save the data on your smartphone:
- weight (👤)
 - BMI
 - body fat percentage (👤)
 - the percentage of water (💧)
 - the percentage of muscle (👤)
 - bone mass in kg (👤)

Attention!
Before taking a new reading, always make sure that the bluetooth connection on your smartphone is active.
If the display shows "E1", that means that the smartphone did not recognize the scale. In this case check:

- that bluetooth is active on the phone
- that the two devices are between 1 and 10 metres apart
- that your smartphone uses a bluetooth 4.0 connection

If the App does not display the data detected by the scale:

- Log out of the App and log back in. After a few minutes your smartphone will show the data detected by your device.

For more information on the LAICA app, visit the website www.laica.com.

TECHNICAL CHARACTERISTICS
Capacity: 180 kg / 397 lb / 28 st 5 lb
Division: 100 g / 0.2 lb
4 - 1.5 V AAA alkaline batteries
Maximum tolerance: +/-1% +1 division
Display LCD reverse
Environmental conditions for use: +10°C +40°C RH 85%
Environmental storage conditions: from +10°C +50°C RH 85%

MANUAL

- Handle the scale with care, it is a precision instrument. Therefore, avoid subjecting the scale to extreme temperature variations, humidity, shocks, dust and direct sunlight.

Instructions and guarantee

- Do not drop or hit the unit.
 - Keep the unit away from heat sources.
 - **Do not jump on the scale.**
 - **Only clean the scale with a soft, dry cloth.**
 - **Never use chemicals or abrasives.**
 - Make sure liquids never penetrate the body of the scale.
- This device is not waterproof. Do not leave or use the scale in wet environments. Humidity greater then 85% or contact with water or other liquids could cause corrosion and compromise proper scale operation.

 Important! Read the instructions carefully

 This product is meant for domestic use.
The conformity, certified by the CE-marking on the device, refers to the EEC Directive 2004/108 on electromagnetic compatibility.
This product is an electronic device that was tested to assure, according to the current technical knowledge, that it does not interfere with other nearby located equipment (electromagnetic compatibility) and that it is safe if used according to the advices given in the instructions for use.
In case of operation anomalies, do not keep using the device and, if necessary, refer to the manufacturer directly.

DISPOSAL PROCEDURE (Dir. 2012/19/Eu-WEEE)

 The symbol on the bottom of the device indicates the separated collection of electric and electronic equipment.
At the end of life of the device, do not dispose it as mixed solid municipal waste,  but dispose it referring to a specific collection centre located in your area or returning it to the distributor, when buying a new device of the same type to be used with the same functions.

If the appliance to be disposed of is less than 25 cm, it can be returned to a retail location that is over 400 m² without having to purchase a new, similar device.
This procedure of separated collection of electric and electronic devices is carried out forecasting a European environmental policy aiming at safeguarding, protecting and improving environment quality, as well as avoiding potential effects on human health due to the presence of hazardous substances in such equipment or to an improper use of the same or of parts of the same.
Caution! The wrong disposal of electric and electronic equipment may involve sanctions.

WARRANTY
This equipment is covered by a 2-year warranty starting from the purchase date, which has to be proven by means **of the stamp or signature of the dealer and by the fiscal payment receipt that shall be kept enclosed to this document.** This warranty period complies with the laws in force and applies only in case the consumer is an individual. Laica products are designed for domestic use only; therefore, their use for commercial purposes is forbidden.
The warranty covers only production defects and does not apply, if the damage is caused by accidental events, misuse, negligence or improper use of the product. Use the supplied accessories only; the use of other accessories invalidate the warranty cover.
Do not open the equipment for any reason whatsoever; the opening or the tampering of the equipment, definitely invalidate the warranty.
Warranty does not apply to parts subject to tear and wear along the use, as well as to batteries, when they are included.
Warranty decays after 2 years from the purchase date. In this case, technical service interventions will be carried out against payment. Information on technical service interventions, under warranty terms or against payment, may be requested by writing at info@laica.com. Repair or replacement of products covered by the warranty terms are free of charge. In case of faults, refer to the dealer; DO NOT send anything directly to LAICA.
All interventions covered by warranty (including product replacement or the replacement of any of its parts) shall not be meant to extend the duration of the original warranty period of the replaced product.
The manufacturer cannot be made liable for any damages that may be caused, directly or indirectly, to people, things and pets resulting from the inobservance of the provisions contained in the suitable user manual, especially with reference to the instructions related to installation, use and maintenance of the equipment. Being constantly involved in improving its own products, Laica reserves the right to modify its product, totally or partially, with reference to any production need, without that this involves any responsibility whatsoever by Laica or its dealers.

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